

PAJUNK®

Manometer



For CSF Pressure Measurement



Spinal Manometer

As the manufacturer of the world's first atraumatic spinal needle, Pajunk takes the diagnostic spinal tap as an integral part of its core competences. With the Sprotte Lumbar needle and our Manometer, Pajunk provides the complete equipment for diagnostic CSF sampling and pressure measurement.

The Pajunk spinal manometer allows for measurements up to 54 cm H₂O (with the additional extension tubing)




1. Solid, transparent tubing

→ Allows for easy measurements up to 34cm H₂O



2. Easy assembly

→ Easy and secure assembly
→ Available in Luer and NRFit



3. Extension tubing

→ Allowing for measurements up to 54cm H₂O
→ Can be purchased separately

Spinal Manometer

Description	Item no. (Luer)	Item no. (NRFit)	PU
3 - 34cm H ₂ O	001151-38G	001163-38F	10

Manometer Extension Tube

Description	Item no.	PU
20cm H ₂ O	001152-38G	10

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

Do you have any questions? Our customer service & sales support will be glad to assist you.

Products distributed in Australia and New Zealand by

SURIMEX 15 Redland Drive T (613) 8833 3804
Mitcham VIC 3132 E sales@surimex.com.au
Australia

www.surimex.com.au